PLATELIA™ CMV IgG

96 TESTS

72680

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE AND QUANTITATIVE DETERMINATION OF
IgG-CLASS ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM

TABLE OF CONTENTS

1. INTENDED USE
2. SUMMARY AND EXPLANATION OF THE TEST
3. PRINCIPLE OF THE TEST
4. KIT CONTENTS AND REAGENT PREPARATION
5. STORAGE AND STABILITY OF REAGENTS
6. PRECAUTIONS
7. TYPE AND STORAGE OF SAMPLE
8. TEST PROCEDURE
9. SCHEME OF TEST PROCEDURE
10. TEST VALIDATION
11. INTERPRETATION OF RESULTS
12. LIMITATIONS OF THE PROCEDURE
13. ANALYTICAL SPECIFICITY
14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY
15. PRECISION
16. "TROUBLE SHOOTING"
17. REFERENCES

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE AND QUANTITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO CYTOMEGALOVIRUS IN HUMAN SERUM

2. SUMMARY AND EXPLANATION OF THE TEST

Cytomegalovirus is a herpes virus transmitted by close human contact. No symptoms of infection are apparent in the majority of cases. However, the virus is very dangerous and may be fatal in immuno-depressed patients. Serum-negative female patients who become infected during pregnancy may transmit the disease to the fetus. In 95% of cases this occurs without symptoms, but some neonates may present jaundice, hepato-splenomegaly and retarded psycho-motorial development. For this reason it is of great importance to determine the immune state of the patient and check for serum conversion. A significant increase in the anti-Cytomegalovirus IgG titer is indicative of a recent infection or reactivation of a latent infection.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay) (1-7).

The antigen, composed of partially purified and inactivated Cytomegalovirus, is bound to the solid phase (8-well strips). The specific immunoglobulins are bound to the antigen through incubation with dilute human serum.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate, composed of human IgG monoclonal antibodies conjugated to horseradish peroxidase.

The unbound conjugate is eliminated, and the peroxidase substrate added.

The blue colour which develops is proportional to the concentration of specific antibodies present in the serum sample. When the enzymatic reaction is interrupted by the addition of a sulphuric acid solution, the yellow colour which develops can be easily read using a microplate reader.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

Bring to room temperature before use.

MT PLATE

MICROPLATE. 12 x 8 wells coated with Cytomegalovirus.

Use: open the package on the opposite side from the code (C plus lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CAL

CALIBRATORS. 5 x 1.6 mL

Contents: Diluted human serum, containing known concentrations of anti-CMV IgG antibodies, in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution. The values can be expressed in IU/ml or EU/ml (see section 11 for conversion table). Calibrator 2 (1 IU/ml or 10 EU/mL) corresponds to the cut-off value and can be used in the qualitative test. The calibrators have the following values: 0.5, 1, 5, 10, 20 IU/mL, obtained by titration against the “**Proposed International Standard WHO**”.

Colour: the colour of the calibrators is proportional to the relative antibody titer.

CONJ

CONJUGATE. 1 x 16 mL.

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

CONTROL IgG - *IgG NEGATIVE CONTROL (PF93910) 1 x 1.6 mL*

INTERCHANGEABLE BETWEEN LOTS

Contents: Human serum not containing anti-CMV IgG antibodies, diluted in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution (Calibrator 0).

Stability: The product is stable up to the expiry date when stored unopened at 2/8°C.

WASH BUF 10x WASH BUFFER 10X (PF93603). 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5%.

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 50x DILUENT 50X (PF93601). 1 x 4.5 mL. For dilution of serum samples.

INTERCHANGEABLE BETWEEN LOTS

Contents: Proteic solution concentrated 50 times, with added phenol 0.05% and Bronidox 0.02%.

Preparation: Dilute the required volume 1:50 in the washing buffer to obtain the diluent ready for use.

SUBS TMB SUBSTRATE (PF93619). 12 mL. Ready for use. **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602). 1 x 16 mL. **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2)

POLYTHENE BAG (1)

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader (wave length 450 or 450/620 nm and 405 nm, with linearity up to OD >= 2000)
- Microplate washer (preferable) able to dispense volumes ranging between 225-375 µl
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	6 weeks at 2/8°C, polythene bag
Calibrators	6 weeks at 2/8°C
Conjugate	6 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	ready for use, 2 weeks at 2/8°C
Rinsing Buffer	2 weeks at 2/8°C, 5 days at 15/30°C
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) the Rinsing Buffer contains detergents
 - b) the conjugate contains phenol
 - c) the substrate is acid
 - d) the controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides; dilute with large amounts of water to eliminate.If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilised after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralised acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30-minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with absorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimised by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyser. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Defrosted samples must be shaken carefully before use. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation.

Human plasma cannot be used in the test.

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Rinsing Buffer 10x (100 mL + 900 mL H₂O).
- Prepare the required diluent for samples by adding 1 part of Diluent 50x to 49 parts of the diluted washing buffer (example: 2 mL + 98 mL of diluted washing buffer).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED calibrators (if possible, in duplicate) in a strip (100 µL in each well). Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µl), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The absorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min. Take the reading again at 405 nm if there are OD higher than 2,000.

9. Scheme of Test Procedure for CYTOMEGALOVIRUS IgG
--

Manual Technique

STEP 1 Place 100 µL of diluted samples/calibrators to the wells of the strips.

Incubate for 45 min. at 37°C

Wash 4 times (300 µl)

STEP 2 Add 100 µL of conjugate to each well

Incubate for 45 min. at 37°C

Wash 4 times (300 µl)

STEP 3 Add 100 µL of Substrate to each well

Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

Read absorbance at 450 nm within 30 min

10. TEST VALIDATION

1. **CONTROL** - (Calibrator 0): the O.D. of the negative control must be < 0.6 times the O.D. of Calibrator 2.
2. The O.D. of Calibrator 2 must be ≥ 0.2 at 450 nm; ≥ 0.16 at 450/620 nm.
3. The O.D. of Calibrator 5 must be over 1.2.

11. INTERPRETATION OF RESULTS

QUANTITATIVE RESULTS

Report the OD of the calibrators on a graph after subtracting the OD of the blank (≤ 0.150). The corresponding titer of the test sample can be obtained by extrapolation.

NOTE: A standard curve must be performed for each run.

If the O.D. of any sample or calibrator is over 2.0, take the reading at 405 nm and multiply the value by 3.

Conversion of the O.D. into units/mL

The anti-Cytomegalovirus IgG can be expressed in EU/mL (arbitrary units) or in IU/mL (proposed anti CMV IgG WHO standard), by extrapolating the results of the 6 calibrators and comparing the OD of the sample with the curve obtained.

Table 1: Calibrator Values in arbitrary and international units

	EU/mL	IU/mL
Calibrator 0	0	0
Calibrator 1	5	0.5
Calibrator 2	10	1
Calibrator 3	50	5
Calibrator 4	100	10
Calibrator 5	200	20

The degree of immunity can be interpreted as follows:

IMMUNE: when the anti-Cytomegalovirus IgG concentration in the sample is > 12 EU/mL or 1.2 IU/mL

NON IMMUNE: when the anti-Cytomegalovirus concentration is < 8 EU/mL or < 0.8 IU/mL

DOUBTFUL: if the result is between the two values. In this case it is advisable to repeat the test in duplicate.

QUALITATIVE RESULTS

Calculate the ratio between the O.D. value of the sample and that of the Cut-off (Calibrator 2). The sample is considered:

Positive: if the ratio is > 1.2

Doubtful: $> 0.8 - < 1.2$

Negative: if the ratio is < 0.8

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

A serum sample obtained during the acute phase of infection, when only IgM antibodies are present, may be negative by this procedure.

Heat-inactivation of the test serum does not give rise to erroneous results in the determination of anti-CMV IgG, as in the case of anti-CMV IgM.

The Cytomegalovirus IgM level should be determined using the Platelia™ CMV IgM kit. Alternatively, a second serum sample obtained 8-14 days later should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of the case history or other diagnostic procedures.

13. ANALYTICAL SPECIFICITY

23 samples which were negative for CMV, but containing IgG antibodies against viruses such as Rubella, Epstein Barr, Herpes Simplex, Measles, Mumps, were tested. In no case did the presence of these antibodies influence the test.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical trial performed in a hospital laboratory, 264 samples were analysed, 47 of which proved negative and 217 positive. The samples were analysed with another commercial immunoenzymatic method: 100% agreement was found between the two methods, both in the positive and the negative samples.

The Platelia™ CMV IgG kit offers 100% sensitivity and specificity.

15. PRECISION

Table 2: Within-run Precision performed on 3 different lots

Cut Off n=21	Batch. N.104	Batch. N. 105	Batch. N. 106
O.D.	0.453	0.374	0.433
CV%	4.62	4.73	4.07

Table 3: "Between-run" Precision

Sample	EU/ml			
	I run	II run	III run	CV%
n.1	12.6	17.5	17.6	18
n.2	78.5	87.5	96.3	10

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert paragraph 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessicant must be pale yellow). Repeat test.
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Colour development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
5. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
6. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
7. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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1. G.B. Wisdom: Enzyme-Immunoassay.
2. Clin. Chem. 22: 1243 (1976).
3. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
4. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
5. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
6. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
7. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
8. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
5. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
6. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
7. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).
8. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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17. REFERENCER

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
5. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
6. M. Musiani et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
7. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980).

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1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. H.O. Kangro: Cytomegalovirus serology: does it give the answer? Serodiagnosis and Immunotherapy 1: 91 (1987).
3. Grint P.C.A. et al.: Screening tests for antibodies to cytomegalovirus: an evaluation of five commercial products. J. Clin. Pathol. 38: 1059 (1985).
4. M. Musiani et al.: Rapid detection of antibodies against cytomegalovirus induced immediate early and early antigens by an enzyme linked immunosorbent assay. J. Clin. Pathol. 37: 122 (1984).
5. F. de Ory et al.: Serological diagnosis of cytomegalovirus infections: comparison of six commercial methods of ELISA. Serodiagnosis and Immunotherapy 2: 423 (1988).
6. P. Dal Monte et al.: Cytomegalovirus umano. Diagnosis 2: 67 (1990).
7. R. Zieglermaier et al.: ELISA. la Ricerca Clin. Lab. 10: 83 (1980)

The other languages which are required in conformity to the European Directive can be obtained from your local Bio-Rad agent.

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